

In-gel Digestion of Excised Protein Spots (Silver Stained Gels) with Destaining

1. Washing/Destaining

- a. Wash gel twice with dd(doubly deionized) H₂O for 15 mins
- b. Excise bands/spots (cut as close to the band/spot as possible to minimize excess Polyacrylamide gel material), and cut into 1 mm cubes and place them in an Eppendorf tube (Polypropylene tube only-- Polypropylene Eppendorf tube, Sarstadt 72.690 & 72.699 (1-800-257-5101)).
- c. Dry samples in a speed vacuum at around 56~60°C for about 45 min, the gel pieces must shrink and be completely dehydrated.
- d. Add 40 µl of 1:1 solution of 30 mM Potassium Ferricyanide and 100 mM Sodium Thiosulfate, incubate at room temperature for 15 min.
- e. Pull off solution and discard, add 40 µl of 100 mM Ambicand incubate at room temperature for 15 min.
- f. Repeat 1e) 3 times
- g. Pull off solution and discard, dry the samples in a speed vacuum as in 1c)

2. Reduction/Alkylation

- a. Remove samples from the speed vacuum and let cool.
- b. Add 40 µl of 10 mM DTT/100 mM Ambic and incubate at 56°C in a water bath or thermocycler for 45 min.
- c. Remove samples and let cool.
- d. Pull off solution and discard, immediately add 40 µl of 55 mM IAA/100 mM Ambic, then incubate at room temperature for 30 min in the dark.
- e. Pull off solution and wash with 40 µl of 100 mM Ambic, then incubate at room temperature for 5 min.
- f. Add 40 µl of acetonitrile to make 1:1 solution of Ambic and ACN, then incubate at room temperature for 15 min.
- g. Pull off solution and discard, dry gel pieces in speed vacuum as in 1c).

3. Digestion/Extraction

- a. Add 40 μ l or enough trypsin solution to cover gel pieces and incubate at 4°C for 45 mins (use ice bath or thermocycler, add more solution if pieces absorb all the liquid).
- b. Pull off excess solution and discard, add 40 μ l (or enough to cover gel pieces) same buffer but without trypsin and incubate at 37°C for 16 hrs (overnight).
- c. Pull off supernatant and save at 3°C.
- d. To gel add 20 μ l of 25 mM Ambic and incubate at room temperature for 15 min.
- e. Add 20 μ l acetonitrile to make 1:1 solution of Ambic/ACN and incubate for 15 min.
- f. Pull off the supernatant and combine it with the one from c).
- g. To gel pieces add 20 μ l of 5% Formic Acid, then incubate at room temperature for 15 min.
- h. Add 20 μ l acetonitrile to make 1:1 solution of ACN/Formic Acid, then incubate for 15 min.
- i. Repeat 3f), 3g), 3h), 3f).
- j. To pooled supernatant add 10 mM DTT to give final concentration of 1 mM DTT .
- k. Completely dry the supernatant, which is the digestion extracts, in speed vacuum.
- l. Resuspend the digestion extracts in 15 μ l~20 μ l of 5% Formic Acid for MS or MS/MS analysis.

Note:

- Larger spots/bands may require more solution. Please adjust the volume accordingly.
- For MALDI analysis, Ziptips(pipette tips with RP-C18 at the tip) should be used to clean up the in-gel digest samples to get rid of salts and detergents prior to spotting the sample on the MALDI state. Routine sensitivities on realworld samples are in the 10~50 fmol rage. Single digit femto-mole level can be reached occasionally. For a routine ziptip protocol, see here.

Some background information:

Mann et. al. have published a paper in Anal Chemistry (1996 vol.68 pg.850-858), where they provide a alternate silver staining method that is MS compatible. However, this staining method is prone to high background, describe here a Rabilloud silver staining method that is better and is also mass spec compatible. Recently Gharahdaghi et al. published a paper in Electrophoresis (1999 vol.20 pg.601-605) that demonstrates a method which increases sensitivity after destaining (remove of silver ions) of the gels pieces prior to enzymatic digestion. So the protocol for enzymatic digestion from silver stained spots is a modified version of Mann's protocol. This protocol has been used with routine real world protein sample analysis from gels in the 10~50ng and higher range.